



FOXP3 Monoclonal Antibody

Product Code	CSB-MA925351
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Purified recombinant Human FOXP3 protein fragments expressed in E.coli.
Raised In	Mouse
Species Reactivity	Human
Specificity	Transfected Only
Tested Applications	ELISA,WB,ICC;Recommended dilution:WB:1:500-1:5000,ICC:1:50-1:200
Relevance	Probable transcription factor.Plays a critical role in the control of immune response.
Form	Ascites containing with 0.02% sodium azide and 50% glycerol.
Isotype	IgG1
Clonality	Monoclonal
Alias	AIID; AIID; DIETER; DIETER; Forkhead box P3; Forkhead box protein P3; foxp3; foxp3; FOXP3_HUMAN; FOXP3delta7; Immune dysregulation polyendocrinopathy enteropathy X linked; Immunodeficiency polyendocrinopathy enteropathy X nked; IPEX; IPEX; JM2; MGC141961
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	7D12
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.