



PTPRC Monoclonal Antibody, Purified

| | |
|----------------------------|--|
| Product Code | CSB-MA135182 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Raised In | Mouse |
| Species Reactivity | Human |
| Tested Applications | ELISA,IF,FC;Recommended dilution:IF:1:100-1:200 |
| Relevance | <p>4AA9 reacts with CD45, a 180-220 kDa leukocyte common antigen (LCA). CD45 antigen is expressed at high levels on all hematopoietic cells including T and B lymphocytes, monocytes, granulocytes, NK cells and dendritic cells, but is not expressed on non-hematopoietic cells. CD45 antibody has also been reported to react weakly with mature blood erythrocytes and platelets. CD45 is a protein tyrosine phosphatase receptor that is critically important for T and B cell antigen receptor-mediated activation.</p> <p>*Target of immunosuppressive antibody treatment.</p> <ol style="list-style-type: none">1. Knapp, W et al., eds. (1989) Leucocyte Typing IV: White Cell Differentiation Antigens, Oxford University Press, New York.2. Trowbridge, IS and Thomas, ML. (1993) Annu. Rev. Immunol. 12:85.3. Pulido, R et al., (1988) J. Immunol. 140:3851.4. Roach, T et al., (1997) Curr. Biol. 7:408. |
| Form | Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide |
| Isotype | IgG2a |
| Clonality | Monoclonal |
| Product Type | Monoclonal Antibody |
| Accession NO. | 4AA9 |