



# CD3E Monoclonal Antibody, FITC/PE/PE-Cy5 Conjugated

<b>Product Code</b>	CSB-MA566015
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Raised In</b>	Mouse
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, IF, FC; Recommended dilution: IF: 1:100-1:200
<b>Form</b>	Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide and 0.2% (w/v) BSA
<b>Isotype</b>	IgG2a/IgG1/IgG2a
<b>Clonality</b>	Monoclonal
<b>Product Type</b>	Monoclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Clone No.</b>	5G4
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.