



CD4/CD8/CD3 Monoclonal Antibody, FITC/PE/PE-Cy5 Conjugated

Product Code	CSB-MA826490
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Raised In	Mouse
Species Reactivity	Human
Tested Applications	FC
Form	Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide and 0.2% (w/v) BSA
Isotype	IgG2b/IgG1/IgG2a
Clonality	Monoclonal
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	7E9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.