



# SELL Monoclonal Antibody, FITC Conjugated

<b>Product Code</b>	CSB-MA289436
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Raised In</b>	Mouse
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, IF, FC; Recommended dilution: IF: 1:100-1:200
<b>Relevance</b>	<p>4AE56 reacts with L-selectin, an 80 kDa leukocyte-endothelial cell adhesion molecule 1 (LECAM-1). CD62L is expressed on most peripheral blood B cells, T cells, some NK cells, monocytes and granulocytes. CD62L mediates lymphocyte homing to high endothelial venules of peripheral lymphoid tissue and leukocyte rolling on activated endothelium at inflammatory sites.</p> <ol style="list-style-type: none"><li>1. Bevilacqua, M et al ( 1991 ) : Cell 67:233.</li><li>2. Polley, MJ et al ( 1991 ) , Proc. Natl. Acad. Sci. USA 88:6224.</li><li>3. Smith, CW et al ( 1991 ) J Clin. Invest . 87:609.</li><li>4. Picker, LJ et al ( 1991 ) Cell 66:921.</li></ol>
<b>Form</b>	Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide and 0.2% (w/v) BSA
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Product Type</b>	Monoclonal Antibody