



# SELP Monoclonal Antibody, PE Conjugated

<b>Product Code</b>	CSB-MA442328
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Raised In</b>	Mouse
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, FC; Recommended dilution: FC: 10 µl/Test
<b>Relevance</b>	<p>?4AW12 reacts with P-selectin , a platelet activation dependent granule-external membrane protein (PADGEM). CD62P is expressed on platelets, megakaryocytes and endothelial cell surface and is upregulated on activated platelets. This molecule mediates rolling of platelets on endothelial cells and rolling of leukocytes on the surface of activated endothelial cells.</p> <p>?</p> <p>1. Schlossman, SL et al., eds. (1995) Leucocyte Typing V : White Cell Differentiation Antigens, Oxford University Press, New York. eds.</p>
<b>Form</b>	Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide and 0.2% (w/v) BSA
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Product Type</b>	Monoclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)