



# IgG1 Monoclonal Antibody;FITC Conjugated

<b>Product Code</b>	CSB-MA255382
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Raised In</b>	Mouse
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	FC
<b>Relevance</b>	This antibody can be used as a mouse IgG1 isotype control in flow cytometry and other applications.
<b>Form</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Product Type</b>	Monoclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Clone No.</b>	22H9
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.