



IgG2a Monoclonal Antibody;FITC Conjugated

Product Code	CSB-MA081360
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Raised In	Mouse
Species Reactivity	Human
Tested Applications	ELISA,IF,FC;Recommended dilution:IF:1:100-1:200
Relevance	This antibody can be used as a mouse IgG2a isotype control in flow cytometry and other applications.
Form	Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide and 0.2% (w/v) BSA
Isotype	IgG2a
Clonality	Monoclonal
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	17C11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.