



# IgG2a Monoclonal Antibody;PE Conjugated

<b>Product Code</b>	CSB-MA072925
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Raised In</b>	Mouse
<b>Tested Applications</b>	ELISA,IF,FC;Recommended dilution:IF:1:100-1:200
<b>Relevance</b>	This antibody can be used as a mouse IgG2a isotype control in flow cytometry and other applications.
<b>Form</b>	Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide and 0.2% (w/v) BSA
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>Product Type</b>	Monoclonal Antibody
<b>Clone No.</b>	3D1
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.