



IgG2b Monoclonal Antibody;FITC Conjugated

Product Code	CSB-MA993519
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Raised In	Mouse
Species Reactivity	Human
Tested Applications	ELISA,IF,FC;Recommended dilution:IF:1:100-1:200
Relevance	This antibody can be used as a mouse IgG2b isotype control in flow cytometry and other applications.
Form	0.2% BSA and 0.1% Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Isotype	IgG2b
Clonality	Monoclonal
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	28E9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.