



IgG2b Monoclonal Antibody;PE Conjugated

Product Code	CSB-MA982879
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Raised In	Mouse
Tested Applications	ELISA,IF,FC;Recommended dilution:IF:1:100-1:200
Relevance	This antibody can be used as a mouse IgG2b isotype control in flow cytometry and other applications.
Form	Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide and 0.2% (w/v) BSA
Isotype	IgG2b
Clonality	Monoclonal
Product Type	Monoclonal Antibody
Clone No.	1E1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.