



CD3E Monoclonal Antibody, PE Conjugated

| | |
|----------------------------|---|
| Product Code | CSB-MA634374 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Raised In | Mouse |
| Species Reactivity | Human |
| Tested Applications | ELISA,IF,FC;Recommended dilution:IF:1:100-1:200 |
| Relevance | 4AT3 reacts with CD3, a 20-26 kDa molecule, which is expressed on all mature T lymphocytes (approximately 60-80% of normal human peripheral blood lymphocytes), NK-T cells and some thymocytes. CD3 associated with the T-cell receptor α/β or γ/δ dimer also plays a role in T-cell activation and signal transduction during antigen recognition.* CD3 mAb is used to count mature T cells in blood. 1. Ledbetter, J et al. (1981) J Exp Med 153:310-323.2. Haynes, B et al. (1986) Leukocyte Typing II: Human T lymphocytes. New York: Sprimger-Verlag; 3-30.3. Clevers H et al. (1988). Ann Rev Immunol 6 : 629-662.4. Barclay, N et al. (1993). The Leucocyte Factsbook. Academic Press. San Diego. |
| Form | Phosphate-buffered solution, pH 7.4, containing 0.09% sodium azide and 0.2% (w/v) BSA |
| Isotype | IgG2a |
| Clonality | Monoclonal |
| Product Type | Monoclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Clone No. | 28C2 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |