



# RPA2 Monoclonal Antibody

<b>Product Code</b>	CSB-MA282022
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Immunogen</b>	Purified recombinant Human RPA32/RPA2 protein fragments expressed in E.coli.
<b>Raised In</b>	Mouse
<b>Species Reactivity</b>	Human
<b>Specificity</b>	This antibody detects endogenous levels of RPA32/RPA2 and does not cross-react with related proteins.
<b>Tested Applications</b>	ELISA,WB,ICC,IP;Recommended dilution:WB:1:500-1:5000,ICC:1:50-1:200IP:1:200-1:2000
<b>Form</b>	Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.
<b>Purification Method</b>	Affinity purified
<b>Isotype</b>	IgG2b
<b>Clonality</b>	Monoclonal
<b>Alias</b>	60S acidic ribosomal protein P1; AA409079; AI325195; AU020965; HSSB; ik:tdsubc_2g1; M(2)21C; MGC137236; OTTHUMP00000004008; p32; p34; RCJMB04_6d17 replication protein A2, 32kDa; REPA 2; REPA1; REPA2; Replication factor A protein 2;
<b>Product Type</b>	Monoclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Clone No.</b>	10H10
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.