



SWAP70 Monoclonal Antibody

Product Code	CSB-MA948113
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Purified recombinant Human SWAP70 protein fragments expressed in E.coli.
Raised In	Mouse
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB;Recommended dilution:WB:1:500-1:5000
Relevance	Phosphatidylinositol 3,4,5-trisphosphate-dependent guanine nucleotide exchange factor (GEF) which,independently of RAS, transduces signals from tyrosine kinase receptors to RAC. It also mediates signaling of membrane ruffling.
Form	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol
Purification Method	Affinity purified
Isotype	IgG2b
Clonality	Monoclonal
Alias	HSPC321, SWAP-70, FLJ39540, KIAA0640
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	11B7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.