



AKT1 Monoclonal Antibody

Product Code	CSB-MA011037
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant fragment of Human AKT1 (C-terminus)
Raised In	Mouse
Species Reactivity	Human,Mouse,Pig,Rat
Tested Applications	ELISA,WB;Recommended dilution:WB:1:500-1:5000
Form	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Purification Method	Antigen Affinity Purified
Isotype	IgG2b
Clonality	Monoclonal
Alias	RAC-alpha serine/threonine-protein kinase PKB alpha RAC-PK-alpha proto-oncogene c-Akt protein kinase B alpha rac protein kinase alpha
Product Type	Monoclonal Antibody