



# ALB Monoclonal Antibody

<b>Product Code</b>	CSB-MA076752
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Immunogen</b>	Recombinant fragment, corresponding to amino acids of Human Albumin
<b>Raised In</b>	Mouse
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA, WB, IHC; Recommended dilution: WB: 1:500-1:5000, IHC: 1:50-1:500
<b>Form</b>	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG2b
<b>Clonality</b>	Monoclonal
<b>Alias</b>	serum albumin albumin (32 AA) albumin (AA 34) growth-inhibiting protein 20 cell growth inhibiting protein 42
<b>Product Type</b>	Monoclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Clone No.</b>	10C5
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.