



ACLY Monoclonal Antibody

Product Code	CSB-MA177768
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant fragment, corresponding to amino acids 734aa-983aa of Human ACLY (C-terminus).
Raised In	Mouse
Species Reactivity	Human,Mouse,Rat,Pig
Tested Applications	ELISA,WB,ICC,FC;Recommended dilution:WB:1:500-1:5000,ICC:1:50-1:200
Form	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Purification Method	Antigen Affinity Purified
Isotype	IgG2a
Clonality	Monoclonal
Alias	ATP-citrate synthase ATP citrate synthase citrate cleavage enzyme ATP-citrate (pro-S-)-lyase
Product Type	Monoclonal Antibody