



CARM1 Monoclonal Antibody

Product Code	CSB-MA256061
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant fragment, corresponding to amino acids 351aa-609aa of Human CARM1(C-terminus).
Raised In	Mouse
Species Reactivity	Human,Rat,Mouse,Pig
Tested Applications	ELISA,WB;Recommended dilution:WB:1:500-1:5000
Form	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Purification Method	Antigen Affinity Purified
Isotype	IgG1
Clonality	Monoclonal
Alias	histone-arginine methyltransferase CARM1 protein arginine N-methyltransferase 4
Product Type	Monoclonal Antibody