



CBX5 Monoclonal Antibody

Product Code	CSB-MA994557
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant full length protein (191aa) of Human CBX5 .
Raised In	Mouse
Species Reactivity	Human,Mouse,Pig,Rat,Dog
Tested Applications	ELISA,WB,ICC;Recommended dilution:WB:1:500-1:5000,ICC:1:50-1:200
Form	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Purification Method	Antigen Affinity Purified
Isotype	IgG2b
Clonality	Monoclonal
Alias	chromobox protein homolog 5 HP1-ALPHA HP1Hs alpha antigen p25 HP1 alpha homolog heterochromatin protein 1-alpha heterochromatin protein 1 homolog alpha chromobox homolog 5 (HP1 alpha homolog, Drosophila)
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	27E2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.