



CHD3 Monoclonal Antibody

Product Code	CSB-MA936397
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant fragment, corresponding to amino acids 193aa-479aa of Human CHD3(C-terminus).
Raised In	Mouse
Species Reactivity	Human,Pig,Dog,Mouse,Rat
Tested Applications	ELISA,WB,IHC,ICC,FC;Recommended dilution:WB:1:500-1:5000,IHC:1:50-1:500,ICC:1:50-1:200
Form	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Purification Method	Antigen Affinity Purified
Isotype	IgG1
Clonality	Monoclonal
Alias	chromodomain-helicase-DNA-binding protein 3 hZFH CHD-3 zinc finger helicase ATP-dependent helicase CHD3 mi-2 autoantigen 240 kDa protein zinc-finger helicase (Snf2-like)
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	5B6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.