



CRTC1 Monoclonal Antibody

Product Code	CSB-MA049428
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant fragment, corresponding to amino acids 25-228aa of Human CRTC1(N-terminus).
Raised In	Mouse
Species Reactivity	Human,Pig,Mouse,Rat
Tested Applications	ELISA,WB,ICC,FC;Recommended dilution:WB:1:500-1:5000,ICC:1:50-1:200
Form	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Purification Method	Antigen Affinity Purified
Isotype	IgG2b
Clonality	Monoclonal
Alias	CREB-regulated transcription coactivator 1 TORC-1 transducer of CREB protein 1 mucoepidermoid carcinoma translocated 1 mucoepidermoid carcinoma translocated protein 1 transducer of regulated cAMP response element-binding protein 1
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	29H11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.