



DsbA Monoclonal Antibody

Product Code	CSB-MA808670
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Full Length (Escherichia coli) Disulfide Oxidoreductase, isolated from E. coli.
Raised In	Mouse
Species Reactivity	E.coli
Tested Applications	ELISA, WB; Recommended dilution: WB: 1:500-1:5000
Form	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol
Purification Method	Antigen Affinity Purified
Clonality	Monoclonal
Alias	hypothetical protein
Product Type	Monoclonal Antibody
Immunogen Species	Escherichia coli