



EPCAM Monoclonal Antibody

Product Code	CSB-MA908159
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant fragment of Human EPCAM (N-terminus).
Raised In	Mouse
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB: 1:500-1:5000
Form	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol
Purification Method	Antigen Affinity Purified
Isotype	IgG1
Clonality	Monoclonal
Alias	epithelial cell adhesion molecule epithelial glycoprotein 314 human epithelial glycoprotein-2 cell surface glycoprotein Trop-1 adenocarcinoma-associated antigen tumor-associated calcium signal transducer 1 major gastrointestinal
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	28G6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.