



FAS Monoclonal Antibody

| | |
|----------------------------|---|
| Product Code | CSB-MA266022 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Immunogen | Recombinant fragment, corresponding to amino acids 211aa-329aa of Human CD95(C-terminus). |
| Raised In | Mouse |
| Species Reactivity | Human |
| Tested Applications | ELISA,WB;Recommended dilution:WB:1:500-1:5000 |
| Form | Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Alias | tumor necrosis factor receptor superfamily member 6 Fas AMA FAS 827dupA CD95 antigen FASLG receptor apoptosis antigen 1 Delta Fas/APO-1/CD95 APO-1 cell surface antigen apoptosis-mediating surface antigen FAS tumor nec |
| Product Type | Monoclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |