



FAS Monoclonal Antibody

Product Code	CSB-MA266022
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant fragment, corresponding to amino acids 211aa-329aa of Human CD95(C-terminus).
Raised In	Mouse
Species Reactivity	Human
Tested Applications	ELISA,WB;Recommended dilution:WB:1:500-1:5000
Form	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Purification Method	Antigen Affinity Purified
Isotype	IgG1
Clonality	Monoclonal
Alias	tumor necrosis factor receptor superfamily member 6 Fas AMA FAS 827dupA CD95 antigen FASLG receptor apoptosis antigen 1 Delta Fas/APO-1/CD95 APO-1 cell surface antigen apoptosis-mediating surface antigen FAS tumor nec
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	1H5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.