



FOXO1 Monoclonal Antibody

Product Code	CSB-MA075284
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant fragment, corresponding to amino acids 420-649aa of Human FOXO1(C-terminus).
Raised In	Mouse
Species Reactivity	Human,Mouse
Tested Applications	ELISA,WB,ICC;Recommended dilution:WB:1:500-1:5000,ICC:1:50-1:200
Form	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Purification Method	Antigen Affinity Purified
Isotype	IgG2a
Clonality	Monoclonal
Alias	forkhead box protein O1 forkhead box protein O1A forkhead, Drosophila, homolog of, in rhabdomyosarcoma
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	30B10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.