



GAPDH Monoclonal Antibody

Product Code	CSB-MA072309
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant fragment, corresponding to amino acids 72aa-248aa of Human GAPDH (Internal).
Raised In	Mouse
Species Reactivity	Human,Pig,Rat,Mouse
Tested Applications	ELISA,WB;Recommended dilution:WB:1:500-1:5000
Form	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Purification Method	Antigen Affinity Purified
Isotype	IgG1
Clonality	Monoclonal
Alias	glyceraldehyde-3-phosphate dehydrogenase aging-associated gene 9 protein peptidyl-cysteine S-nitrosylase GAPDH glyceraldehyde 3-phosphate dehydrogenase
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)