



GAPDH Monoclonal Antibody

| | |
|----------------------------|---|
| Product Code | CSB-MA072309 |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Immunogen | Recombinant fragment, corresponding to amino acids 72aa-248aa of Human GAPDH (Internal). |
| Raised In | Mouse |
| Species Reactivity | Human,Pig,Rat,Mouse |
| Tested Applications | ELISA,WB;Recommended dilution:WB:1:500-1:5000 |
| Form | Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Alias | glyceraldehyde-3-phosphate dehydrogenase aging-associated gene 9 protein peptidyl-cysteine S-nitrosylase GAPDH glyceraldehyde 3-phosphate dehydrogenase |
| Product Type | Monoclonal Antibody |
| Immunogen Species | Homo sapiens (Human) |
| Clone No. | 12H5 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |