



HK2 Monoclonal Antibody

Product Code	CSB-MA105781
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant fragment, corresponding to amino acids 270-435aa of Human HK2(N-terminus).
Raised In	Mouse
Species Reactivity	Human,Pig,Mouse,Rat
Tested Applications	ELISA,WB,FC;Recommended dilution:WB:1:500-1:5000
Form	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Purification Method	Antigen Affinity Purified
Isotype	IgG1
Clonality	Monoclonal
Alias	hexokinase-2 HK II hexokinase type II hexokinase-2, muscle muscle form hexokinase mi-2 autoantigen 240 kDa protein zinc-finger helicase (Snf2-like)
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	8E8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.