



HSPA4 Monoclonal Antibody

Product Code	CSB-MA559162
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant fragment, corresponding to amino acids 96-289aa of Human HSPA1A
Raised In	Mouse
Species Reactivity	Human,Dog,Pig,Mouse,Rat
Tested Applications	ELISA,WB;Recommended dilution:WB:1:500-1:5000
Form	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Purification Method	Antigen Affinity Purified
Isotype	IgG1
Clonality	Monoclonal
Alias	heat shock 70 kDa protein 1A/1B HSP70-1/HSP70-2 HSP70.1/HSP70.2 heat shock 70kD protein 1A heat shock-induced protein heat shock 70 kDa protein 1/2 dnaK-type molecular chaperone HSP70-1
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	10C12
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.