



# HSPA4 Monoclonal Antibody

<b>Product Code</b>	CSB-MA559162
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Immunogen</b>	Recombinant fragment, corresponding to amino acids 96-289aa of Human HSPA1A
<b>Raised In</b>	Mouse
<b>Species Reactivity</b>	Human,Dog,Pig,Mouse,Rat
<b>Tested Applications</b>	ELISA,WB;Recommended dilution:WB:1:500-1:5000
<b>Form</b>	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Alias</b>	heat shock 70 kDa protein 1A/1B HSP70-1/HSP70-2 HSP70.1/HSP70.2 heat shock 70kD protein 1A heat shock-induced protein heat shock 70 kDa protein 1/2 dnaK-type molecular chaperone HSP70-1
<b>Product Type</b>	Monoclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)