



IKZF1 Monoclonal Antibody

Product Code	CSB-MA782104
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant fragment, corresponding to amino acids 249-415aa of Human IKAROS(C-terminus).
Raised In	Mouse
Species Reactivity	Human,Mouse,Dog
Tested Applications	ELISA,WB;Recommended dilution:WB:1:500-1:5000
Form	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Purification Method	Antigen Affinity Purified
Isotype	IgG1
Clonality	Monoclonal
Alias	DNA-binding protein Ikaros CLL-associated antigen KW-6 Ikaros (zinc finger protein) ikaros family zinc finger protein 1 lymphoid transcription factor LyF-1 zinc finger protein, subfamily 1A, 1 (Ikaros
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	12E1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.