



NFKBIA Monoclonal Antibody

Product Code	CSB-MA448062
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant fragment, corresponding to amino acids of Human NFKBIA
Raised In	Mouse
Species Reactivity	Human,Pig,Mouse,Rat
Tested Applications	ELISA,WB;Recommended dilution:WB:1:500-1:5000
Form	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Purification Method	Antigen Affinity Purified
Clonality	Monoclonal
Alias	NF-kappa-B inhibitor alpha ikB-alpha IkappaBalpha kappa-B-alpha nuclear factor of kappa light chain gene enhancer in B-cells major histocompatibility complex enhancer-binding protein MAD3
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	9F9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.