



JMJD6 Monoclonal Antibody

Product Code	CSB-MA184549
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant fragment, corresponding to amino acids 11-248aa of Human JMJD6(N-terminus).
Raised In	Mouse
Species Reactivity	Human,Dog,Mouse,Rat
Tested Applications	ELISA,WB,ICC;Recommended dilution:WB:1:500-1:5000,ICC:1:50-1:200
Form	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Purification Method	Antigen Affinity Purified
Isotype	IgG1
Clonality	Monoclonal
Alias	bifunctional arginine demethylase and lysyl-hydroxylase JMJD6 lysyl-hydroxylase JMJD6 phosphatidylserine receptor jmjC domain-containing protein 6 histone arginine demethylase JMJD6 peptide-lysine 5-dioxygenase JMJD6 jumonji doma
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	5G6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.