



# KDM1A Monoclonal Antibody

<b>Product Code</b>	CSB-MA581953
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Immunogen</b>	Recombinant fragment (Human)
<b>Raised In</b>	Mouse
<b>Species Reactivity</b>	Human,Mouse,Rat,Pig
<b>Tested Applications</b>	ELISA,WB;Recommended dilution:WB:1:500-1:5000
<b>Form</b>	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
<b>Purification Method</b>	Antigen Affinity Purified
<b>Clonality</b>	Monoclonal
<b>Alias</b>	lysine-specific histone demethylase 1A lysine (K)-specific demethylase 1 BRA35-HDAC complex protein BHC110 lysine-specific histone demethylase 1 amine oxidase (flavin containing) domain 2 FAD-binding protein BRA35-HDAC complex, 110
<b>Product Type</b>	Monoclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)