



KRT7 Monoclonal Antibody

Product Code	CSB-MA031700
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant fragment, , corresponding to amino acids 256aa-415aa of Human KRT7 (C-terminus).
Raised In	Mouse
Species Reactivity	Human,Dog,Mouse,Rat
Tested Applications	ELISA,WB,ICC;Recommended dilution:WB:1:500-1:5000,ICC:1:50-1:200
Form	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Purification Method	Antigen Affinity Purified
Isotype	IgG2b
Clonality	Monoclonal
Alias	keratin, type II cytoskeletal 7 CK-7 keratin-7 sarcolectin cytokeratin 7 cytokeratin-7 type-II keratin Kb7 type II mesothelial keratin K7 keratin, 55K type II cytoskeletal keratin, simple epithelial type I, K7
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	5A3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.