



CCAR2 Monoclonal Antibody

Product Code	CSB-MA264953
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant fragment, corresponding to amino acids 1aa-265aa of Human KIAA1967(N-terminus).
Raised In	Mouse
Species Reactivity	Human,Mouse,Rat,Pig,Dog
Tested Applications	ELISA,WB;Recommended dilution:WB:1:500-1:5000
Form	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Purification Method	Antigen Affinity Purified
Isotype	IgG1
Clonality	Monoclonal
Alias	p30 DBC protein DBC.1 deleted in breast cancer 1 deleted in breast cancer gene 1 protein
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	29B3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.