



TP53 Monoclonal Antibody

Product Code	CSB-MA189580
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant fragment of Human TP53 (C-terminus)
Raised In	Mouse
Species Reactivity	Human,Dog,Mouse,Rat
Tested Applications	ELISA,WB;Recommended dilution:WB:1:500-1:5000
Form	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Purification Method	Antigen Affinity Purified
Isotype	IgG1
Clonality	Monoclonal
Alias	cellular tumor antigen p53 antigen NY-CO-13 phosphoprotein p53 p53 tumor suppressor transformation-related protein 53
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)