



POU5F1 Monoclonal Antibody

Product Code	CSB-MA586746
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant full length protein (352aa) of Mouse POU5F1
Raised In	Mouse
Species Reactivity	Human
Tested Applications	ELISA,WB,ICC,FC;Recommended dilution:WB:1:500-1:5000,ICC:1:50-1:200
Form	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Purification Method	Antigen Affinity Purified
Isotype	IgG1
Clonality	Monoclonal
Alias	POU domain, class 5, transcription factor 1 octamer-binding protein 3 octamer-binding protein 4 POU domain transcription factor OCT4 octamer-binding transcription factor 3 octamer-binding transcription factor-3 POU-type homeodoma
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	13C12
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.