



PPP1CA Monoclonal Antibody

Product Code	CSB-MA590830
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant fragment, corresponding to amino acids 2-286aa of Human PPP1CA (N-terminus).
Raised In	Mouse
Species Reactivity	Human,Pig,Rat,Mouse,Rabbit
Tested Applications	ELISA,WB,ICC;Recommended dilution:WB:1:500-1:5000,ICC:1:50-1:200
Form	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Purification Method	Antigen Affinity Purified
Isotype	IgG1
Clonality	Monoclonal
Alias	serine/threonine-protein phosphatase PP1-alpha catalytic subunit protein phosphatase 1, catalytic subunit, alpha isoform serine/threonine protein phosphatase PP1-alpha 1 catalytic subunit
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	10D4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.