



PPARG Monoclonal Antibody

Product Code	CSB-MA930213
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant fragment, corresponding to amino acids 254-478aa of Human PPARG(C-terminus).
Raised In	Mouse
Species Reactivity	Human,Dog,Pig,Rat,Mouse
Tested Applications	ELISA,WB,ICC;Recommended dilution:WB:1:500-1:5000,ICC:1:50-1:200
Form	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Purification Method	Antigen Affinity Purified
Isotype	IgG1
Clonality	Monoclonal
Alias	peroxisome proliferator-activated receptor gamma PPAR gamma PPAR-gamma nuclear receptor subfamily 1 group C member 3 peroxisome proliferator-activated receptor gamma 1 peroxisome proliferator-activated nuclear receptor gamma variant
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	23F2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.