



PRMT6 Monoclonal Antibody

Product Code	CSB-MA912675
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant full length protein (375aa) of Human PRMT6.
Raised In	Mouse
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC,ICC;Recommended dilution:WB:1:500-1:5000,IHC:1:50-1:500,ICC:1:50-1:200
Form	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Purification Method	Antigen Affinity Purified
Isotype	IgG1
Clonality	Monoclonal
Alias	protein arginine N-methyltransferase 6 HMT1 hnRNP methyltransferase-like 6 histone-arginine N-methyltransferase PRMT6 heterogeneous nuclear ribonucleoprotein methyltransferase-like protein 6
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	8B11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.