



PPP4C Monoclonal Antibody

Product Code	CSB-MA994759
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant full length protein (307aa) of Human PPP4C.
Raised In	Mouse
Species Reactivity	Human,Dog,Mouse,Bovine,Rat
Tested Applications	ELISA,WB,IHC,ICC;Recommended dilution:WB:1:500-1:5000,IHC:1:50-1:500,ICC:1:50-1:200
Form	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Purification Method	Antigen Affinity Purified
Isotype	IgG2a
Clonality	Monoclonal
Alias	serine/threonine-protein phosphatase 4 catalytic subunit PP-X protein phosphatase X, catalytic subunit protein phosphatase 4 (formerly X), catalytic subunit
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	4B7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.