



# PPP4C Monoclonal Antibody

<b>Product Code</b>	CSB-MA994759
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Immunogen</b>	Recombinant full length protein (307aa) of Human PPP4C.
<b>Raised In</b>	Mouse
<b>Species Reactivity</b>	Human,Dog,Mouse,Bovine,Rat
<b>Tested Applications</b>	ELISA,WB,IHC,ICC;Recommended dilution:WB:1:500-1:5000,IHC:1:50-1:500,ICC:1:50-1:200
<b>Form</b>	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>Alias</b>	serine/threonine-protein phosphatase 4 catalytic subunit PP-X protein phosphatase X, catalytic subunit protein phosphatase 4 (formerly X), catalytic subunit
<b>Product Type</b>	Monoclonal Antibody
<b>Immunogen Species</b>	Homo sapiens (Human)