



PDHA1 Monoclonal Antibody

Product Code	CSB-MA275281
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant fragment , corresponding to amino acids 131aa-389aa of Human PDHA1 (C-terminus).
Raised In	Mouse
Species Reactivity	Human,Dog,Rat,Mouse,Pig
Tested Applications	ELISA,WB,ICC;Recommended dilution:WB:1:500-1:5000,ICC:1:50-1:200
Form	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Purification Method	Antigen Affinity Purified
Isotype	IgG1
Clonality	Monoclonal
Alias	pyruvate dehydrogenase E1 component subunit alpha, somatic form, mitochondrial PDHE1-A type I pyruvate dehydrogenase complex, E1-alpha polypeptide 1
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	12C9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.