



SMAD5 Monoclonal Antibody

Product Code	CSB-MA052999
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant fragment, corresponding to amino acids 239aa-372aa of Human SMAD5 (C-terminus).
Raised In	Mouse
Species Reactivity	Human,Pig,Rat,Mouse
Tested Applications	ELISA,WB,ICC,FC;Recommended dilution:WB:1:500-1:5000,ICC:1:50-1:200
Form	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Purification Method	Antigen Affinity Purified
Isotype	IgG1
Clonality	Monoclonal
Alias	mothers against decapentaplegic homolog 5 SMA- and MAD-related protein 5 SMAD, mothers against DPP homolog 5 MAD, mothers against decapentaplegic homolog 5 mothers against decapentaplegic, drosophila, homolog of, 5
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)