



SMC1A Monoclonal Antibody

Product Code	CSB-MA064269
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant fragment, corresponding to amino acids 2aa-351aa of Human SMC1A(N-terminus).
Raised In	Mouse
Species Reactivity	Human,Pig,Mouse,Rat
Tested Applications	ELISA,WB,IHC,FC;Recommended dilution:WB:1:500-1:5000,IHC:1:50-1:500
Form	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Purification Method	Antigen Affinity Purified
Isotype	IgG1
Clonality	Monoclonal
Alias	structural maintenance of chromosomes protein 1A SMC-1A SMC-1-alpha SMC protein 1A segregation of mitotic chromosomes 1 SMC1 structural maintenance of chromosomes 1-like 1 SMC1 (structural maintenance of chromosomes 1, yeast)-lik
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)