



SOX2 Monoclonal Antibody

Product Code	CSB-MA046266
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant full length protein (319aa) of Mouse Sox2.
Raised In	Mouse
Species Reactivity	Human,Pig,Mouse,Rat
Tested Applications	ELISA,WB;Recommended dilution:WB:1:500-1:5000
Form	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Purification Method	Antigen Affinity Purified
Isotype	IgG2b
Clonality	Monoclonal
Alias	transcription factor SOX-2 transcription factor SOX2 SRY-related HMG-box gene 2
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)