



TAB1 Monoclonal Antibody

Product Code	CSB-MA964181
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant fragment, corresponding to amino acids of Human TAB1
Raised In	Mouse
Species Reactivity	Human,Dog,Mouse,Pig,Rat
Tested Applications	ELISA,WB;Recommended dilution:WB:1:500-1:5000
Form	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Purification Method	Antigen Affinity Purified
Clonality	Monoclonal
Alias	TGF-beta-activated kinase 1 and MAP3K7-binding protein 1 TAK1-binding protein 1 transforming growth factor beta-activated kinase-binding protein 1 mitogen-activated protein kinase kinase kinase 7-interacting protein 1
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	28C4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.