



TBL1XR1 Monoclonal Antibody

Product Code	CSB-MA953191
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Immunogen	Recombinant full length protein (514aa) of Human TBLR1.
Raised In	Mouse
Species Reactivity	Human,Mouse,Rat,Pig,Dog
Tested Applications	ELISA,WB,IHC,ICC;Recommended dilution:WB:1:500-1:5000,IHC:1:50-1:500,ICC:1:50-1:200
Form	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50%,glycerol
Purification Method	Antigen Affinity Purified
Isotype	IgG1
Clonality	Monoclonal
Alias	F-box-like/WD repeat-containing protein TBL1XR1 TBL1-related protein 1 transducin beta-like 1X-related protein 1 nuclear receptor co-repressor/HDAC3 complex subunit nuclear receptor co-repressor/HDAC3 complex subunit TBLR1
Product Type	Monoclonal Antibody
Immunogen Species	Homo sapiens (Human)
Clone No.	3H4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.